




NORMAL LABORATORY VALUES, Page 1 of 2

IMPORTANT NOTES: Each commercial laboratory has its own set of "normal" values, called "Normal Range" or "Reference Range" on your lab report. These values depend on the equipment or method used. Compare your results to the range shown on your lab report. Results that are "out of range" may not represent a problem. Your test results can be affected by several factors, including your age or gender, if you are pregnant, the time of day when the sample was taken, active infections, stage of HIV disease, and food (some test samples need to be taken after you have fasted - not eaten anything - for several hours). Where normal values for men and women are different, they are indicated as W for women and M for men. **Discuss "out of range" results with your health care provider.** The table below compares the units used in the United States with the "Système International d'Unités (SI units), a metric system used in many parts of the world. The last column, "To Convert US to SI Units," is the factor to multiply US lab values to convert them to SI units. To convert SI units to US units, divide the SI value by the conversion factor. See page 2 for a terminology list.

Laboratory Test	Normal Range in US Units	Normal Range in SI Units	To Convert US to SI Units
ALT (Alanine aminotransferase)	W 7-30 units/liter M 10-55 units/liter	W 0.12-0.50 µkat/liter M 0.17-0.92 µkat/liter	x 0.01667
Albumin	3.1 - 4.3 g/dl	31 - 43 g/liter	x 10
Alkaline Phosphatase	W 30-100 units/liter, M 45-115 units/liter	W 0.5-1.67 µkat/liter, M 0.75-1.92 µkat/liter	x 0.01667
Amylase (Serum)	53-123 units/liter	0.88-2.05 nkat/liter	x 0.01667
AST (Aspartate aminotransferase)	W 9-25 units/liter M 10-40 units/liter	W 0.15-0.42 µkat/liter M 0.17-0.67 µkat/liter	x 0.01667
Basophils	0-3% of lymphocytes	0.0-0.03 fraction of white blood cells	x 0.01
Bilirubin - Direct	0.0-0.4 mg/dl	0-7 µmol/liter	x 17.1
Bilirubin - Total	0.0-1.0 mg/dl	0-17 µmol/liter	
Blood pressure	Normal: 120/70 to 120/80 millimeters of mercury (mmHg). Top number is systolic pressure, when heart is pumping. Bottom number is diastolic pressure when heart is at rest. Blood pressure can be too low (hypotension) or too high (hypertension).		No conversion
C peptide	0.5-2.0 ng/ml	0.17-0.66 nmol/liter	x 0.33
Calcium, serum	8.5-10.5 mg/dl	2.1-2.6 mmol/liter	x 0.25
Calcium, urine	0-300 mg/24hr	0.0-7.5 mmol/24hr	x 0.025
CO ₂ (Bicarbonate)	20 - 32 mmol/L	20 - 32 mmol/L	No conversion
Chloride	95 - 108 mmol/L	95 - 108 mmol/L	No conversion
Cholesterol, Total & LDL Desirable Marginal High Very High HDL: Desirable Moderate Low (heart risk)	Total Cholesterol <200 mg/dL	LDL cholesterol <100 mg/dL	x0.02586
	200-239 mg/dL	100 - 159 mg/dL	
	>239 mg/dL	160 - 189 mg/dL	
	--	>190 mg/dL	
	>60 mg/dL	>1.55 mmol/liter	
	40 - 60 mg/dL	1.03 - 1.55 mmol/liter	
	<40 mg/dL	<1.03 mmol/liter	
Cortisol: serum free (urine)	0-25 µg/dl (depends on time of day)	0-690 nmol/liter	x 27.59
	20-70 µg/dl	55-193 nmol/24hr	x 2.759
Creatine kinase	W 40-150 units/liter, M 60-400 units/liter	W 0.67-2.50 µkat/liter, M 1.00-6.67 µkat/liter	x 0.01667
Creatinine (urine)	W 0.6 - 1.8 g/day M 0.8 - 2.4 g/day	W 5.3-15.9 mmol/day M 7.1-21.2 mmol/day	x 88.4
DHEA	W 130-980 ng/dl M 180-1250 ng/dl	W 4.5-34.0 nmol/liter M 6.24-43.3 nmol/liter	x 0.03467
DHEA Sulfate	W Pre-menopause: 12-535 µg/dl W Post-menopause: 30-260 µg/dl M 10-619 µg/dl	W Pre-menopause: 120-5350 µg/liter W Post-menopause: 300-2600 µg/liter M 100-6190 µg/liter	x 10
Eosinophils	0-8% of white blood cells	0.0-0.8 fraction of white blood cells	x 0.01
Erythrocyte sedimentation rate (Sed Rate)	W ≤ 30 mm/h M ≤ 20 mm/h	W ≤ 30 mm/h M ≤ 20 mm/h	No conversion
Folate	3.1-17.5 ng/ml	7.0-39.7 nmol/liter	x 2.266
Glucose, urine	<0.05 g/dl	<0.003 mmol/liter	x 0.05551
Glucose, plasma	70-110 mg/dl	3.9-6.1 mmol/liter	
Gamma glutamyl transferase (GGT)	W ≤ 45 U/L M ≤ 65 U/L	W ≤ 45 U/L M ≤ 65 U/L	No conversion



NORMAL LABORATORY VALUES, Page 2

Laboratory Test	Normal Range in US Units	Normal Range in SI Units	To Convert US to SI Units
Hematocrit	W 36.0% - 46.0% of red blood cells M 37.0% - 49.0% of red blood cells	W 0.36-0.46 fraction of red blood cells M 0.37-0.49 fraction of red blood cells	x 0.01
Hemoglobin	W 12.0-16.0 g/dl, M 13.0-18.0 g/dl	W 7.4-9.9 mmol/liter, M 8.1-11.2 mmol/liter	x 0.6206
Lactate dehydrogenase (LDH) (total)	≤ 270 U/L	≤ 4.5 µkat/liter	X 0.016667
Lactic acid	0.5-2.2 mmol/liter	0.5-2.2 mmol/liter	No conversion
Leukocytes (WBC)	4.5-11.0x10 ³ /mm ³	4.5-11.0x10 ⁹ /liter	No conversion
Lymphocytes	16-46% of white blood cells	0.16-0.46 fraction of white blood cells	x 0.01
Mean corpuscular hemoglobin (MCH)	25.0-35.0 pg/cell	25.0-35.0 pg/cell	No conversion
Mean corpuscular hemoglobin concentration (MCHC)	31.0-37.0 g/dl	310-370 g/liter	x 10
Mean corpuscular volume (MCV)	W 78-102 µm ³ M 78-100 µm ³	W 78-102 fl M 78-100 fl	No conversion
Monocytes	4-11% of white blood cells	0.04-0.11 fraction of white blood cells	x 0.01
Neutrophils	45-75% of white blood cells	0.45-0.75 fraction of white blood cells	x 0.01
Phosphorus	2.5 – 4.5 mg/dL	0.81 – 1.45 mmol/L	X 0.323
Platelets (Thrombocytes)	130 – 400 x 10 ³ /µL	130 – 400 x 10 ⁹ /L	No conversion
Potassium	3.4-5.0 mmol/liter	3.4-5.0 mmol/liter	No conversion
Red Blood Cell Count (RBC)	W 3.9–5.2x10 ⁶ /µL M 4.4–5.8x10 ⁶ /µL	W 3.9–5.2 x 10 ¹² /L M W 4.4 – 5.8 x 10 ¹² /L	No conversion
Sodium	135-145 mmol/liter	135-145 mmol/liter	No conversion
Testosterone, total (morning sample)	W 6-86 ng/dl M 270-1070 ng/dl	W 0.21-2.98 nmol/liter M 9.36-37.10 nmol/liter	x 0.03467
Testosterone, Unbound	Age 20-40 W 0.6-3.1, M 15.0-40.0 pg/ml Age 41-60 W 0.4-2.5, M 13.0-35.0 pg/ml Age 61-80 W 0.2-2.0, M 12.0-28.0 pg/ml	W 20.8-107.5, M 520-1387 pmol/liter W 13.9-86.7, M 451-1213 pmol/liter W 6.9-69.3, M 416-971 pmol/liter	x 34.67
Triglycerides (fasting)	Normal 40-150 mg/dl Borderline 150-200 mg/dl High 200-500 mg/dl Very High >500 mg/dl	0.45-1.69 mmol/liter 1.69 - 2.26 mmol/liter 2.26 - 5.65 mmol/liter >5.65 mmol/liter	x 0.01129
Urea, plasma (BUN)	8-25 mg/dl	2.9-8.9 mmol/liter	x 0.357
Urinalysis: pH	5.0-9.0	5.0-9.0	No conversion
Specific gravity	1.001-1.035	1.001-1.035	No conversion
WBC (White blood cells, Leukocytes)	4.5-11.0x10 ³ /mm ³	4.5-11.0x10 ⁹ /liter	No conversion

TERMINOLOGY: UNITS:

gram: common measurement of weight. Used in this table: pg (picograms), g (grams), mg (milligrams), etc. per liter

katal (kat): a unit of catalytic activity, used especially in the chemistry of enzymes. Used in this table: µkat (microkatal), nkat (nanokatal) per liter

micrometer (µm): a unit of length. Mean Corpuscular Volume is expressed in cubic micrometers

mole: also "gram molecular weight," a quantity based on the atomic weight of the substance. Many test results in the Système Internationale are expressed as the number of moles per liter. In US units, these measurements are usually in grams per liter. Used in this table: mmol (millimoles), µmol, (micromoles), nmol (nanomoles), pmol (picomoles) per liter

Some units of measurement include the following fractions and multipliers:

mega (M): 10 ⁶ or x1,000,000	milli (m): 10 ⁻³ or ÷1,000
kilo (k): 10 ³ or x1,000	micro (µ): 10 ⁻⁶ or ÷1,000,000
deca or deka: 10 ¹ or x10	nano (n): 10 ⁻⁹ or ÷1,000,000,000
deci (d): 10 ⁻¹ or ÷10	pico (p): 10 ⁻¹² or ÷1,000,000,000,000

Reviewed September 30, 2014